

ANALYTICAL SPECIFICATION

Issued by: K. Bäckström

Edition no: 0403AG

Approved by: A.N.de Belder

Supersedes: 0403AF

Valid from: 12-2016

Dextran sulfate 5 kDa HS

Catalogue No. DB004

Description: A polyanionic derivative of dextran with a weight average molecular weight of approximately 5000. Supplied as the sodium salt as a white to off white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-6000	02030
Number average molecular weight (Mn)	to be noted	02030
Sulfur content	16.0-20.0 %	02011
pH	5.0-7.5	02009
Free sulfates	≤ 2%	02013
Loss on drying	≤ 10%	02018
