

ANALYTICAL SPECIFICATION

Issued by: G. Jandikova

Edition no: 0203AE

Approved by: Shams Ghanem

Supersedes: 0203AD

Valid from: 02-2020

FITC-dextran 10 kDa
Catalogue No. FD10

Description: A fluorescent derivative of a dextran fraction with an average molecular weight of approximately 10000. It is supplied as a yellow powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow to dark orange	02037
Appearance, form	Powder	02037
Weight average molecular weight (M_w)	9000-11000	02030
Number average molecular weight (M_n)	To be noted	02030
Degree of substitution	0.001-0.008	02002
Loss on drying	≤ 10%	02018
Free FITC	Pass	02005
