

ANALYTICAL SPECIFICATION

Issued by: G. Jandikova

Edition no: 0207AE

Approved by: Shams Ghanem

Supersedes: 0207AD

Valid from: 02-2020

FITC-dextran 150 kDa

Catalogue No. FD150

Description: A fluorescent derivative of a dextran fraction with an average molecular weight of approximately 150000. It is supplied as a yellow to dark orange powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow or light orange to dark orange or yellow-orange or orange-yellow	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	140000-160000	02030
Number average molecular weight (Mn)	To be noted	02030
Loss on drying	≤ 10%	02018
Degree of substitution	0.001-0.008	02002
Free FITC	Pass	02005
