

ANALYTICAL SPECIFICATION

Issued by: E. Hallgren

Edition no: 0202AD

Approved by: C. de Belder Tesseus

Supersedes: 0202AC

Valid from: 03-2018

FITC-dextran 4 kDa

Catalogue No. FD4

Description: A fluorescent derivative of a dextran fraction with an average molecular weight of approximately 4000. It is supplied as a yellow powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow to dark orange	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-5000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.008	02002
Loss on drying	≤ 10%	02018
Free FITC	Pass	02005
