

ANALYTICAL SPECIFICATION

Issued by:	Gabriela Jandikova	Edition no: 0405AE
Approved by:	Charlotte de Belder Tesséus	Supersedes: 0405AD
		Valid from: 10-2022

Dextran sulfate 10 HS

Catalogue No. DB008

Description: A polyanionic derivative of dextran with a weight average molecular weight of approximately 10000. Supplied as the sodium salt as a white to off white powder readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	9000-16000	02030
Number average molecular weight (Mn)	To be noted	02030
Loss on drying	≤ 10%	02018
Sulphur content	16-20 %	02050
рН	5.0-7.5	02009
Free sulfates	≤ 2%	02013