

ANALYTICAL SPECIFICATION

Issued by: K. Bäckström **Edition no:** 0408AD

Approved by: A.N.de Belder **Supersedes:** 0408AC

Valid from: 12-2016

Dextran sulfate 20 LS

Catalogue No. DB013

Description: A polyanionic derivative of dextran with a weight average molecular weight of

approximately 20000. Supplied as the sodium salt as a white to off white powder

which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	16000-30000	02030
Number average molecular weight (Mn)	To be noted	02030
Loss on drying	≤ 10%	02018
Sulphur content	8-13 %	02011
Specific optical rotation	+90° to +140°	02008
рН	5.0-7.5	02009
Free sulfates	≤ 2%	02013