

ANALYTICAL SPECIFICATION

Tentative

Issued by: G. Jandikova **Edition no:** 2505AB

Approved by: Charlotte de B Tesséus **Supersedes:** 2505AA

Valid from: 04-2020

Lysine-dextran 40

Catalogue No. LD40

Description: A Lysine derivative of a dextran with a weight average molecular weight of

approximately 40000. Supplied as a white to off white powder which is readily

soluble in water.

| TEST/CHARACTERISTIC | LIMITS | TEST METHOD |
|---|--------------------|-------------|
| Appearance, colour | White to off white | 02037 |
| Appearance, form | Powder | 02037 |
| Weight average molecular weight (Mw) | 30000-50000 | 02030 |
| Number average molecular weight (Mn) | To be noted | 02030 |
| Degree of substitution | > 0.001 | 02041 |
| Specific optical rotation | To be noted | 02008 |
| рН | 5.0-8.0 | 02009 |
| Free lysine | Pass | 02042 |
| Loss on drying | ≤ 10% | 02018 |