

## ANALYTICAL SPECIFICATION

## **Tentative**

**Issued by:** G. Jandikova **Edition no:** 2508AB

**Approved by:** Charlotte de B Tesséus **Supersedes:** 2504AA

**Valid from:** 04-2020

## Lysine-dextran 500

Catalogue No. LD500

**Description:** A Lysine derivative of a dextran with a weight average molecular weight of

approximately 500000. Supplied as a white to off white powder which is readily

soluble in water.

| TEST/CHARACTERISTIC                     | LIMITS             | TEST METHOD |
|---|--------------------|-------------|
| Appearance, colour                      | White to off white | 02037       |
| Appearance, form                        | Powder             | 02037       |
| Weight average<br>molecular weight (Mw) | 400000-700000      | 02030       |
| Number average<br>molecular weight (Mn) | To be noted        | 02030       |
| Degree of substitution                  | > 0.001            | 02041       |
| Specific optical rotation               | To be noted        | 02008       |
| рН                                      | 5.0-8.0            | 02009       |
| Free lysine                             | Pass               | 02042       |
| Loss on drying                          | ≤ 10%              | 02018       |