

ANALYTICAL SPECIFICATION

Tentative

Issued by: G. Jandikova **Edition no:** 2507AA

Approved by: Charlotte de B Tesséus Supersedes: -

Valid from: 10-2020

Lysine-dextran 150

Catalogue No. LD150

Description: A Lysine derivative of a dextran with a weight average molecular weight of

approximately 150000. Supplied as a white to off white powder which is readily

soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	To be noted	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	> 0.001	02041
Specific optical rotation	To be noted	02008
рН	5.0-8.0	02009
Free lysine	Pass	02042
Loss on drying	≤ 10%	02018