

## ANALYTICAL SPECIFICATION

**Tentative** 

Issued by: Gabriela Jandikova Edition no: 1502AD

**Approved by:** Charlotte de Belder Tesséus **Supersedes:** 1502AC

Valid from: 04-2021

## **DEAE-dextran 20**

Catalogue No. DD20

Description: A diethylaminoethyl-dextran fraction with an average molecular weight of

approx. 20000. Supplied as the sodium salt as a white powder which is readily

soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	15000-30000	02030
Number average molecular weight (Mn)	To be noted	02030
Nitrogen content	2.5-4.5 %	02023
рН	5.0-6.5	02009
Loss on drying	≤10%	02018