

ANALYTICAL SPECIFICATION

Issued by:	Shams Ghanem	Edition no: 1501AC
Approved by:	A.N.de Belder	Supersedes: 1501AB

Valid from: 07-2017

DEAE-dextran 4

Catalogue No. DD4

Description: A diethylaminoethyl-dextran fraction with an average molecular weight of approx. 4000. Supplied as the sodium salt as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-6000	02030
Number average molecular weight (Mn)	To be noted	02030
Nitrogen content	2-5 %	02023
рН	5.0-6.5	02009
Loss on drying	≤ 10%	02018