

ANALYTICAL SPECIFICATION

Issued by: Shams Ghanem **Edition no:** 1503AC

Approved by:A.N. de BelderSupersedes: 1503AB

Valid from: 07-2017

DEAE-dextran 40

Catalogue No. DD40

Description: A diethylaminoethyl-dextran fraction with an average molecular weight of

approx. 40000. Supplied as the sodium salt as a white powder which is readily

soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	30000-55000	02030
Number average molecular weight (Mn)	To be noted	02030
Nitrogen content	2-5 %	02023
рН	5.0-6.5	02009
Loss on drying	≤10%	02018