

ANALYTICAL SPECIFICATION

Issued by: Gabriela Jandikova Edition no: 0607AD

Approved by: C. de Belder Tesséus **Supersedes:** 0607AC

Valid from: 09-2020

Blue dextran 20

Catalogue No. BD20

Description: Blue dextran 20 is prepared from dextran with an average molecular weight of

approximately 20000 Dalton by reacting with Cibacron blue F3GA. It is supplied

as a blue powder, soluble in water.

CAS No. -

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Blue to dark blue	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	15000-25000	02030
Number average molecular weight (Mn)	To be noted	02030
Free Cibacron blue	Pass	02007
Degree of substitution (mmol/g)	0.01-0.08	02004
Loss on drying	≤ 10%	02018