

ANALYTICAL SPECIFICATION

Issued by:	Gabriela Jandikova		Edition no: 0602AE
Approved by:	Charlotte de Belder Tesséus		Supersedes: 0602AD
			Valid from: 05-2022
Blue dextran 5 Catalogue No. BD5			
Description:	A derivative of a dextran fraction with an average molecular weight of approximately 5000. It is supplied as a blue powder which is slowly soluble in water.		
CAS No.	87915-38-6		
TEST/CHARACTERISTIC		LIMITS	TEST METHOD
Appearance, colour		Blue to dark blue	02037
Appearance, form		Powder	02037
Weight average molecular weight (Mw)		4000-6000	02030
Number average molecular weight (Mn)		To be noted	02030
Degree of substitution (mmol/g)		0.01-0.08	02004
Free Cibacron blue		Pass	02007
Loss on drying		≤ 10%	02018