

## ANALYTICAL SPECIFICATION

**Tentative** 

**Issued by:** G. Jandikova **Edition no:** 1403AD

**Approved by:** C. de Belder Tesséus **Supersedes:** 1403AC

**Valid from:** 06-2022

## FITC-dextran sulfate 10

Catalogue No. FDSS10

**Description:** A FITC derivative of dextran sulfate 10 (DSS). Supplied as a yellow to pale yellow

powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	7500-15000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.020	02002
Sulphur content	16-20 %	02011
Free FITC	≤ 100 ppm	02005
Loss on drying	≤ 10%	02018