

ANALYTICAL SPECIFICATION

Tentative

Issued by: G. Jandikova **Edition no:** 1404AD

Approved by:C. de Belder TesséusSupersedes: 1404AC

Valid from: 06-2022

FITC-dextran sulfate 4

Catalogue No. FDSS4

Description: A FITC derivative of dextran sulfate 4 (DSS). Supplied as a yellow to pale yellow

powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-6000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.020	02002
Sulphur content	16-20 %	02011
Free FITC	≤ 100 ppm	02005
Loss on drying	≤ 10%	02018