

ANALYTICAL SPECIFICATION

Issued by: G. Jandikova Edition no: 1401AE

Approved by: C. de Belder Tesséus **Supersedes:** 1401AD

Valid from: 06-2022

FITC-dextran sulfate 40

Catalogue No. FDSS40

Description: A FITC derivative of dextran sulfate 40 (DSS). Supplied as a yellow to pale

yellow powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	30000-50000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.020	02002
Sulphur content	16-19 %	02011
Free FITC	≤ 100 ppm	02005
Loss on drying	≤ 10%	02018