

ANALYTICAL SPECIFICATION

Tentative

Issued by: G. Jandikova **Edition no:** 2902AE

Approved by: Charlotte de B Tesséus **Supersedes:** 2902AD

Valid from: 06-2022

FITC-lysine-dextran 4

Catalogue No. FLD4

Description: A fluorescent derivative of lysine-dextran with an average molecular weight of

approx. 4000 which is supplied as a yellow powder readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow to dark orange	02037
Appearance, form	Powder or flakes	02037
Weight average molecular weight (Mw)	2000-6000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.03	02002
Free FITC	≤ 100 ppm	02005
Loss on drying	≤ 10%	02018