

ANALYTICAL SPECIFICATION

Issued by: G. Jandikova **Edition no:** 0205AI

Approved by: Charlotte de Belder Tesséus **Supersedes:** 0205AH

Valid from: 06-2022

FITC-dextran 40

Catalogue No. FD40

Description: A fluorescent derivative of a dextran fraction with an average molecular weight

of approximately 40000. It is supplied as a yellow powder which is readily soluble

in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Light yellow to dark yellow to dark orange	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	35000-45000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.020	02002
Free FITC	≤ 100 ppm	02005
Loss on drying	≤ 10%	02018