

ANALYTICAL SPECIFICATION

Issued by: Gabriela Jandikova Edition no: 0508AD

Approved by:Raffaello PapadakisSupersedes: 0508AC

Valid from: 08-2022

TRITC-dextran 2000

Catalogue No. TD2000

Description: A fluorescent derivative of a dextran fraction with an average molecular weight

of approximately 2000000. It is supplied as a pink to red powder which is

readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Pink to purple	02037
Appearance, form	Powder or crystals	02037
Weight average molecular weight (Mw)	Made from dextran 2000	02030
Degree of substitution	0.001-0.008	02003
Free TRITC	≤ 100 ppm	02006
Loss on drying	≤ 10 %	02018