

ANALYTICAL SPECIFICATION

Issued by: Gabriela Jandikova Edition no: 0502AF

Approved by:Charlotte de Belder TesséusSupersedes: 0502AE

Valid from: 08-2022

TRITC-dextran 4

Catalogue No. TD4

Description: A fluorescent derivative of a dextran fraction with an average molecular weight

of approximately 4000. It is supplied as a red powder which is readily soluble in

water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Pink	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-5000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.008	02003
Free TRITC	≤ 100 ppm	02006
Loss on drying	≤ 10 %	02018