

ANALYTICAL SPECIFICATION

Issued by: Gabriela Jandikova Edition no: 3404AB

Approved by: Charlotte de B Tesséus **Supersedes:** 3404AA

Valid from: 08-2022

TRITC-lysine-dextran 4

Catalogue No. TLD4

Description: A fluorescent derivative of a Lysine-dextran fraction with an average molecular

weight of approximately 4000 which is supplied as a pink powder readily soluble

in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Pink	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	3000-5000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.03	02003
Free TRITC	≤ 100 ppm	02006
Loss on drying	≤ 10 %	02018