

## ANALYTICAL SPECIFICATION

## Tentative

Issued by:	G. Jandikova	Edition no: 3403AC
Approved by:	Charlotte de B Tesséus	Supersedes: 3403AB
		Valid from: 08-2022

TRITC-lysine-dextran 500 Catalogue No. TLD500

**Description:** A fluorescent derivative of a Lysine-dextran fraction with an average molecular weight of approximately 500000 which is supplied as a pink powder readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	Pink	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	400000-800000	02030
Number average molecular weight (Mn)	To be noted	02030
Degree of substitution	0.001-0.03	02003
Free TRITC	≤ 100 ppm	02006
Loss on drying	≤ 10%	02018