

ANALYTICAL SPECIFICATION

Issued by:	G.Jan
issued by:	G.Jan

Approved by:

idikova

C. de Belder Tesséus

Edition no: 0436AB

Supersedes: 0436AA

Valid from: 05-2023

Dextran sulfate 10 HS, Ph. grade Catalogue No. DS10 HS PG

Description: A polyanionic derivative of dextran with a weight average molecular weight of approximately 10000. Supplied as the sodium salt as a white to off white powder which is readily soluble in water.

TEST/CHARACTERISTIC	2	LIMITS	TEST METHOD
Appearance, colour		White to off white	02037
Appearance, form		Powder	02037
Weight average molecular weight (Mw)		9000-16000	02030
Number average mol	ecular weight (Mn)	To be noted	02030
Sulphur content		16-20 %	02011
рН		5.0-7.5	02009
Free sulfate		≤ 0.2 %	02013
Loss on drying		≤ 7 %	02018
Specific optical rotati	on	+80° to +105°	02008
Turbidity		< 5 NTU	02021
Residual solvents:	EtOH	To be noted	Ph.Eur. 2.4.24
	Formamide	To be noted	Ph.Eur. 2.4.24
Total aerobic microbial count (TAMC)		≤ 10 ³ CFU/g	Ph.Eur. 2.6.12
Total combined yeasts/mould count (TYMC)		≤ 10² CFU/g	Ph.Eur. 2.6.12
Bacterial endotoxins (BET)		< 5 IU/mg	Ph.Eur. 2.6.14

We hereby confirm that no metal catalysts or metal reagents are used in the manufacturing of this product. Therefore, elemental impurities, classified according to ICH Q3D, are unlikely to be present.

We hereby confirm that no class 1 solvent, classified according to Ph.Eur. 5.4 and USP <467> Residual solvents, is used in the manufacturing of this product.