

ANALYTICAL SPECIFICATION

Issued by:

G.Jandikova

Approved by:

C. de Belder Tesséus

Edition no: 0429AC

Supersedes: 0429AB

Valid from: 05-2023

	Valid from: 05-		
		i n sulfate 5 LS, Ph. grade talogue No. DS5 LS PG	
Description:	A polyanionic derivative of dextran with a weight average molecular weight of approximately 5000. Supplied as the sodium salt as a white to off white powder which is readily soluble in water.		
TEST/CHARACTERISTIC		LIMITS	TEST METHOD
Appearance, colour		White to off white	02037
Appearance, form		Powder	02037
Weight average molecular weight (Mw)		3000-6000	02030
Number average molecular weight (Mn)		To be noted	02030
Sulphur content		8.0-13.0 %	02011
рН		5.0-7.5	02009
Free sulfate		≤ 0.2 %	02013
Loss on drying		≤ 7 %	02018
Specific optical rotation		+65° to +95°	02008
Turbidity		< 5 NTU	02021
Residual solve	nts		
EtOH		To be noted (According to ICH Q3C)	Ph.Eur. 2.4.24
Formamide		To be noted (According to ICH Q3C)	Ph.Eur. 2.4.24
Total aerobic microbial count (TAMC)		≤ 10³ CFU/g	Ph.Eur. 2.6.12
Total combined yeasts/mould count (TYMC)		≤ 10² CFU/g	Ph.Eur. 2.6.12
Bacterial endotoxins (BET)		< 5 IU/mg	Ph.Eur. 2.6.14

We hereby confirm that no metal catalysts or metal reagents are used in the manufacturing of this product. Therefore, elemental impurities, classified according to ICH Q3D, are unlikely to be present.

We hereby confirm that no class 1 solvent, classified according to Ph.Eur. 5.4 and USP <467> Residual solvents, is used in the manufacturing of this product.