

## ANALYTICAL SPECIFICATION

**Issued by:** Gabriela Jandikova **Edition no:** 1703AD

Approved by:Charlotte de Belder TesséusSupersedes: 1703AC

**Valid from:** 09-2020

## Polysucrose 1000

Catalogue No. P1000

**Description:** A polysucrose fraction with an average molecular weight of approx. 1000000 Da.

Supplied as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder or crystals or granules	02037
Weight average molecular weight (Mw)	900000-1300000	02030
Number average molecular weight (Mn)	To be noted	02030
Specific optical rotation	+30° to +55°	02008
рН	To be noted	02009
Absorbance	< 0.1	02015
Loss on drying	≤ 10%	02018