

ANALYTICAL SPECIFICATION

Issued by: Gabriela Jandikova Edition no: 1701AD

Approved by: Charlotte de Belder Tesséus **Supersedes:** 1701AC

Valid from: 11-2020

Polysucrose 20

Catalogue No. P20

Description: A polysucrose fraction with an average molecular weight of approx. 20000 Da.

Supplied as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	15000-30000	02030
Number average molecular weight (Mn)	To be noted	02030
Specific optical rotation	To be noted	02008
рН	To be noted	02009
Absorbance	< 0.1	02015
Loss on drying	≤ 10%	02018