

ANALYTICAL SPECIFICATION

Tentative

Issued by: Gabriela Jandikova Edition no: 1705AA

Approved by: Charlotte de Belder Tesséus Supersedes: -

Valid from: 03-2021

Polysucrose 400

Catalogue No.-

Description: A polysucrose fraction with an average molecular weight of approx. 400000 Da.

Supplied as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	300000-500000	02030
Number average molecular weight (Mn)	To be noted	02030
Specific optical rotation	To be noted	02008
рН	To be noted	02009
Absorbance	< 0.2	02015
Loss on drying	≤ 10%	02018