

ANALYTICAL SPECIFICATION

Tentative

Issued by: Gabriela Jandikova

Edition no: 1702AD

Approved by: Charlotte de Belder Tesséus

Supersedes: 1702AC

Valid from: 09-2020

Polysucrose 40
Catalogue No. P40

Description: A polysucrose fraction with an average molecular weight of approx. 40000 Da. Supplied as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	35000-50000	02030
Number average molecular weight (Mn)	To be noted	02030
Specific optical rotation	To be noted	02008
pH	To be noted	02009
Absorbance	< 0.2	02015
Loss on drying	≤ 10%	02018
