

ANALYTICAL SPECIFICATION

Tentative

Issued by: Gabriela Jandikova

Edition no: 1706AA

Approved by: Charlotte de Belder Tesséus

Supersedes: -

Valid from: 08-2023

Polysucrose 50
Catalogue No. P50

Description: A polysucrose fraction with an average molecular weight of approx. 50000 Da. Supplied as a white powder which is readily soluble in water.

TEST/CHARACTERISTIC	LIMITS	TEST METHOD
Appearance, colour	White to off white	02037
Appearance, form	Powder	02037
Weight average molecular weight (Mw)	45000-60000	02030
Number average molecular weight (Mn)	To be noted	02030
Specific optical rotation	To be noted	02008
pH	To be noted	02009
Absorbance	< 0.2	02015
Loss on drying	≤ 10 %	02018
